

WATER BATH LTBAT31-4



WATER BATH LTBAT31-4

LTBAT31-4 WATER BATH

PID temperature controller, large LCD display screen and interface provides for userfriendly operation

Easy to set Adjustable Timer. (1 minute to 5,999 minutes)

Built-in circulating water pump to ensure uniform upper and lower temperature of bath lotion

Maintenance-free operation with easy to clean surface

Both heater and bath chamber are made of corrosion-resistant stainless steel

No angle in bath chamber, easy to clean

The standard stainless-steel bottom plate, helps prevents direct contact by accessories and tubes to heating element and sensors

Drain valve makes emptying of bath water fast and easy for cleaning and moving



SPECIFICATIONS

Model	LTBAT31-4
Capacity	22 L
Power requirements	AC 220 V 50 HZ
Power	1000 W
Temperature Range	RT+5 - 80°C
Temperature fluctuation	±0.2°C
Temperature resolution	0.1°C
Internal dimension	505Wx150Dx330H mm
External dimension	558Wx340Dx342H mm
Timer	1 - 5999 min
Net weight	19 Kg
Porous cover(lid)	6 - hole

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com