

LOW TEMPERATURE CIRCULATOR LTBAT34-16



LOW TEMPERATURE CIRCULATOR LTBAT34-16

LTBAT34-16 LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



SPECIFICATIONS

| | |
|-------------------------|-------------|
| Model | LTBAT34-16 |
| Capacity | 13 L |
| Temperature range | -50-100°C |
| Precision | ±0.2 |
| Interior Dimension | 240x170x200 |
| Electrical Requirements | 220 V 50 Hz |
| Pump (flux) | 8 L/min |
| Power Consumption | 3100 W |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com