

# SHAKING WATER BATH LTBAT34-2



# SHAKING WATER BATH LTBAT34-2

## LTBAT34-2 SHAKING WATER BATH

Microprocessor temperature controller with LCD screen ensures precise and reliable control, easy to operate

Both heater and bath chamber are made of corrosion-resistant stainless steel

No angle in bath chamber, easy to clean

Programming setting function with 7 periods and 9 steps for each period, which means there are 63 programmable steps in total

Easy to set adjustable timer. (1 minute to 5,999 minutes)

Maintenance-free operation with easy to clean surface

Drain valve makes emptying of bath water fast and easy for cleaning and moving



## SPECIFICATIONS

Model	LTBAT34-2
Capacity	47 L
Temperature range	RT+5 - 99°C
Display resolution	0.1°C
Temperature uniformity (at 37°C)	±1.0°C
Shaking speed range	30 - 180 rpm
Amplitude	30mm (Standard) or 40mm (Option)
Interior dimension	618Wx310Dx250H mm
Exterior dimension	823Wx350Dx355H mm
Power requirements	AC 220 V 50 HZ
Power consumption	1650 W



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)