

LOW TEMPERATURE CIRCULATOR LTBAT34-5



LOW TEMPERATURE CIRCULATOR LTBAT34-5

LTBAT34-5 LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	LTBAT34-5
Capacity	22.5 L
Temperature range	RT+5-100°C
Precision	±0.1
Interior Dimension	330x300x150
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	1050 W

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com