

# LOW TEMPERATURE CIRCULATOR LTBAT34-6



# LOW TEMPERATURE CIRCULATOR LTBAT34-6

## LTBAT34-6 LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



### SPECIFICATIONS

Model	LTBAT34-6
Capacity	14.5 L
Temperature range	RT+5-100°C
Precision	±0.2
Interior Dimension	240x170x200
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	1050 W

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)