

LOW TEMPERATURE CIRCULATOR LTBAT36-3



LOW TEMPERATURE CIRCULATOR LTBAT36-3

LTBAT36-3 LOW TEMPERATURE CIRCULATOR

High precision temperature controller



SPECIFICATIONS

Model	LTBAT36-3
Capacity	13 L
Temperature range	-40-100°C
Precision	±0.02
Peristome Dimension	240x170x200
Pump(flux)	15 L/min
Power	2850 W

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com