

# SHAKING WATER BATH LTBAT37-3



# SHAKING WATER BATH LTBAT37-3

## LTBAT37-3 SHAKING WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



## SPECIFICATIONS

Model	LTBAT37-3
Capacity	47 L
Temperature range	RT+5-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1650 W
Interior Dimension	618Wx310Dx250H mm
Exterior Dimension	823Wx350Dx355H mm



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)