

SHAKING WATER BATH LTBAT37-4



SHAKING WATER BATH LTBAT37-4

LTBAT37-4 SHAKING WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



SPECIFICATIONS

| | |
|------------------------|----------------------------------|
| Model | LTBAT37-4 |
| Capacity | 47 L |
| Temperature range | RT+5-100°C |
| Display Resolution | 0.1°C |
| Temperature Uniformity | ±1°C |
| Shaking speed range | 30-150rpm |
| Amplitude | 30mm (Standard) or 40mm (Option) |
| Power Consumption | 1650 W |
| Interior Dimension | 618Wx310Dx250H mm |
| Exterior Dimension | 823Wx350Dx355H mm |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com