

WATER BATH LTBAT38-2



WATER BATH LTBAT38-2

LTBAT38-2 WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



SPECIFICATIONS

Model	LTBAT38-2
Capacity	4 holes $\Phi 92$ mm
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	1000 W
Temperature Stability	$\pm 0.3^{\circ}\text{C}$
Temp Alarm	$\pm 2^{\circ}\text{C}$
Interior dimension	220Wx280Dx150H mm
Exterior dimension	396Wx330Dx260H mm
Timing Range	1~5999 min



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com