

# WATER BATH LTBAT38-3



# WATER BATH LTBAT38-3

## LTBAT38-3 WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



## SPECIFICATIONS

Model	LTBAT38-3
Capacity	6 holes $\Phi$ 92 mm
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	2050 W
Temperature Stability	$\pm$ 0.3°C
Temp Alarm	$\pm$ 2°C
Interior dimension	290Wx490Dx150H mm
Exterior dimension	600Wx390Dx260H mm
Timing Range	1~5999 min

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)