

WATER BATH LTBAT39-1



WATER BATH LTBAT39-1

LTBAT39-1 WATER BATH

Microprocessor control with timing function
Digital display
Audible and visible alarm for over temperature



SPECIFICATIONS

Model	LTBAT39-1
Capacity	11 L
Electrical Requirement	220 V 50 HZ
Power Consumption	500 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	420Wx180Dx150H mm
Exterior dimension	570Wx220Dx260H mm
Timing Range	1~999 min



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com