

WATER BATH LTBAT39-3



WATER BATH LTBAT39-3

LTBAT39-3 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



SPECIFICATIONS

Model	LTBAT39-3
Capacity	12 L
Electrical Requirement	220 V 50 HZ
Power Consumption	400 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	300Wx240Dx160H mm
Exterior dimension	460Wx280Dx190H mm
Timing Range	1~999 min



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com