

# WATER BATH LTBAT39-4



# WATER BATH LTBAT39-4

## LTBAT39-4 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



### SPECIFICATIONS

Model	LTBAT39-4
Capacity	22 L
Electrical Requirement	220 V 50 HZ
Power Consumption	600 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	450Wx300Dx190H mm
Exterior dimension	610Wx340Dx260H mm
Timing Range	1~999 min



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)