

# WATER BATH LTBAT39-5



# WATER BATH LTBAT39-5

## LTBAT39-5 WATER BATH

Microprocessor control with timing function  
Digital display  
Audible and visible alarm for over temperature



### SPECIFICATIONS

Model	LTBAT39-5
Capacity	30 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	600Wx300Dx190H mm
Exterior dimension	760Wx340Dx280H mm
Timing Range	1~999 min

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)