

WATER BATH LTBAT39-6



WATER BATH LTBAT39-6

LTBAT39-6 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



SPECIFICATIONS

Model	LTBAT39-6
Capacity	2.1 L x 3
Electrical Requirement	220 V 50 HZ
Power Consumption	750 W
Temperature Range	RT+5-99°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	150Wx125Dx110H mm
Exterior dimension	490Wx245Dx310H mm
Timing Range	1~999 min



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com