

HIGH TEMPERATURE CIRCULATOR LTBAT40-1



HIGH TEMPERATURE CIRCULATOR LTBAT40-1

LTBAT40-1 HIGH TEMPERATURE CIRCULATOR

Microprocessor control with timing function
Digital display
Audible and visible alarm for over temperature



SPECIFICATIONS

Model	LTBAT40-1
Capacity	22 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-70°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	500Wx300Dx150H mm
Exterior dimension	660Wx340Dx200H mm
Timing Range	1~999 min

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com