

HIGH TEMPERATURE CIRCULATOR LTBAT40-2



HIGH TEMPERATURE CIRCULATOR LTBAT40-2

LTBAT40-2 HIGH TEMPERATURE CIRCULATOR

Microprocessor control with timing function
Digital display
Audible and visible alarm for over temperature



SPECIFICATIONS

Model	LTBAT40-2
Capacity	30 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-70°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	600Wx300Dx190H mm
Exterior dimension	770Wx370Dx280H mm
Timing Range	1~999 min

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com