

# AUTOMATED BLOOD CULTURE SYSTEM LTBCS8-5



# AUTOMATED BLOOD CULTURE SYSTEM

## LTBCS8-5

Used in Clinical laboratory, Medical laboratory.

Also known as Automated Microbial Detection System, Instrumented Blood Culture System, Microbial Detection Instrument Microbial Culture Instrument.

## LTBCS8-5 AUTOMATED BLOOD CULTURE SYSTEM

Four incubator drawers of 60 cell capacity deliver 240 cell capacity

One heating system for each incubator drawer to control temperature

Control module handles up to 64 incubator modules

Support for random loading and unloading, support for incubation of anonymous bottles, positive anonymous and negative anonymous bottles can be identified

Double-barcode design to load/unload bottles simply and avoid operation mistakes

Time to positive (TTP) of more than 90% positive specimens is within 24 hours

Alarm by audio and light flash



## SPECIFICATIONS

Model	LTBCS8-5
Capacity	240
Dimension (LxWxH)	780x650x980 mm
Weight	200 kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)