

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-2



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-2

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB8-2 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	LTCAB8-2
Capacity	12 gallon (45L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	590x460x890 mm
Package Size	710x580x1066 mm
Gross Weight	66 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com