

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-2



# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-2

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## LTCAB8-2 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



### SPECIFICATIONS

Model	LTCAB8-2
Capacity	12 gallon (45L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	590x460x890 mm
Package Size	710x580x1066 mm
Gross Weight	66 Kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)