

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-3



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-3

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB8-3 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	LTCAB8-3
Capacity	22 gallon (83L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	3 pcs
Door	Single / Manual
External Size	600x460x1650 mm
Package Size	720x580x1815 mm
Gross Weight	116 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com