

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-4



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-4

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB8-4 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	LTCAB8-4
Capacity	30 gallon (114L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Double / Manual
External Size	1090x460x1120 mm
Package Size	1210x680x1285 mm
Gross Weight	117 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com