

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-4



# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-4

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.  
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## LTCAB8-4 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



### SPECIFICATIONS

Model	LTCAB8-4
Capacity	30 gallon (114L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Double / Manual
External Size	1090x460x1120 mm
Package Size	1210x680x1285 mm
Gross Weight	117 Kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)