

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-5



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-5

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB8-5 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	LTCAB8-5
Capacity	45 gallon (170L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1650mm
Package Size	1210x580x1815 mm
Gross Weight	155 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com