

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-7



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET LTCAB8-7

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB8-7 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	LTCAB8-7
Capacity	90 gallon (340L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x860x1650 mm
Package Size	1210x980x1815 mm
Gross Weight	243 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com