

STRONG ACID AND ALKALI STORAGE CABINET LTCAB9-2



STRONG ACID AND ALKALI STORAGE CABINET

LTCAB9-2

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB9-2 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali,
Customizable



SPECIFICATIONS

Model	LTCAB9-2
Capacity	12 gallon (45L)
Material	pp
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	590x460x890 mm
Package Size	675x545x1025 mm
Gross Weight	30 Kg

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com