

STRONG ACID AND ALKALI STORAGE CABINET LTCAB9-3




STRONG ACID AND ALKALI STORAGE CABINET

LTCAB9-3

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB9-3 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali,
Customizable



SPECIFICATIONS

Model	LTCAB9-3
Capacity	20 gallon (75L)
Material	pp
Adjustable Shelf	1 pc
Door	Double / Manual
External Size	910x590x890 mm
Package Size	995x675x1035 mm
Gross Weight	46 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com