

STRONG ACID AND ALKALI STORAGE CABINET LTCAB9-5




STRONG ACID AND ALKALI STORAGE CABINET

LTCAB9-5

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB9-5 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali,
Customizable



SPECIFICATIONS

Model	LTCAB9-5
Capacity	30 gallon (114L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1120mm
Package Size	1175x545x1255 mm
Gross Weight	58 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com