

STRONG ACID AND ALKALI STORAGE CABINET LTCAB9-6



STRONG ACID AND ALKALI STORAGE CABINET

LTCAB9-6

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB9-6 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali,
Customizable



SPECIFICATIONS

Model	LTCAB9-6
Capacity	45 gallon (170L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1650mm
Package Size	1175x545x1785 mm
Gross Weight	72 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com