

STRONG ACID AND ALKALI STORAGE CABINET LTCAB9-7



STRONG ACID AND ALKALI STORAGE CABINET

LTCAB9-7

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB9-7 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali,
Customizable



SPECIFICATIONS

Model	LTCAB9-7
Capacity	60 gallon (227L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	860x860x1650 mm
Package Size	945x945x1785 mm
Gross Weight	90 Kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com