

AIR SHOWER BOX LTCBS21-1



AIR SHOWER BOX LTCBS21-1

Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.
Also known as Pass Box.

LTCBS21-1 AIR SHOWER BOX

Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness



SPECIFICATIONS

Model	LTCBS21-1
Clean Air Speed	0.36~0.54m/s
UV Lamp	15Wx1
External Size (WxDxH)	660x570x1000 mm
Internal Size (WxDxH)	500x500x500 mm
Material	304 stainless steel
Interlock	Electronical
HEPA Filter	Efficiency $\geq 99.995\%$ at $0.3\mu\text{m}$
Gross Weight	80 kg
Packing Size(WxDxH)	780x700x1250 mm
Power Supply	220V, 50Hz; 110V, 50/60Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com