

AIR SHOWER BOX LTCBS21-1



AIR SHOWER BOX LTCBS21-1

Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.
Also known as Pass Box.

LTCBS21-1 AIR SHOWER BOX

Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness



SPECIFICATIONS

| | |
|-----------------------|------------------------------|
| Model | LTCBS21-1 |
| Clean Air Speed | 0.36~0.54m/s |
| UV Lamp | 15Wx1 |
| External Size (WxDxH) | 660x570x1000 mm |
| Internal Size (WxDxH) | 500x500x500 mm |
| Material | 304 stainless steel |
| Interlock | Electronical |
| HEPA Filter | Efficiency ≥99.995% at 0.3μm |
| Gross Weight | 80 kg |
| Packing Size(WxDxH) | 780x700x1250 mm |
| Power Supply | 220V, 50Hz; 110V, 50/60Hz |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com