

BENCHTOP CENTRIFUGE LTCEN16-8



BENCHTOP CENTRIFUGE LTCEN16-8

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuges, Laboratory Floor Type Centrifuges, Benchtop Centrifuges, Non Refrigerated High Speed Centrifuges, Laboratory Tabletop Centrifuges.

LTCEN16-8 BENCHTOP CENTRIFUGE

Gentle braking at low speeds

LCD display

Short Spin key for quick centrifugation

Lid is opened automatically after the operation is finished, which avoids overheating of the sample and save processing time

Safety lid-lock



SPECIFICATIONS

Model	LTCEN16-8
Maximum Capacity (No of tubes x Vol.)	6x15ml
Maximum Speed	100-5000
Maximum RCF	3074xg
Time Range	15 S~99 min
Acceleration/ Deceleration	20-80s↑ 20-80s↓
Power	40 W
Power Supply	110V 50/60Hz, 230V-50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
2501615008	Rotor	D12x75 glass tube	6x15ml
2501615009	Rotor	Conical tube with screw cap	8x15 ml
2501615010	Rotor	Conical tube with screw cap	12x5 ml



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com