

MICROCENTRIFUGE LTCEN29-10



MICROCENTRIFUGE LTCEN29-10

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.
Also known as Benchtop Centrifuges, Laboratory Centrifuges.

LTCEN29-10 MICROCENTRIFUGE

Quick release of rotor for easy rotor exchange
Gentle braking system for quick sample processing



SPECIFICATIONS

Model	LTCEN29-10
Maximum Capacity (No of tubes x Vol.)	8x1.5 ml / 2.0 ml
Speed Range	4000 & 6000
Maximum RCF	1200xg & 2000xg
Overall Dimension	135x147x120 mm
Weight (Net/Gross)	0.7 kg
Power Supply	AC100-240V

OPTIONAL ACCESSORIES

Accessory Code	Name	Capacity
2502917008	Fixed-Angle Rotor	8x1.5/2.0ml
2502917009	Strip Rotor	4-8x0.2ml
2502917010	Adapter	0.5-0.6ml tube
2502917011	Adapter	0.2ml tube



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com