

BACTERIAL COLONY COUNTER LTCOU8-2



BACTERIAL COLONY COUNTER LTCOU8-2

LTCOU8-2 BACTERIAL COLONY COUNTER

CMOS integrated circuit and LED digital display.

The special probe counts sensitively and accurately.

It is composed by counter, probe and count pond, etc.

The magnification is variable by adjusting the flectional arm.



SPECIFICATIONS

Model	LTCOU8-2
Phonetic Counting	No
Amplification	5X, 10X
Counter Range	0 ~ 999
Suitable Petri Dishes	φ50-150mm
Character Height	13mm
Lamp Power	16W
External Size(W*D*H)	270*200*70mm
Package Size(W*D*H)	370*270*330mm
Net Weight	1.7kg
Gross Weight	3kg
Power Supply	AC110V, 50/60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com