

# COLONY SPREADER LTCOU8-3



# COLONY SPREADER LTCOU8-3

## LTCOU8-3 COLONY SPREADER

Four speeds options for user's requirement

Aluminum alloy disc can be sterilized by disinfectants and UV

Continuous operation is for daily use

Feasible 70/90mm plate, no need to replace the adapter



## SPECIFICATIONS

Model	LTCOU8-3
Diameter	170 mm
Height	80 mm
Operative Temperature	15°C~30°C
Centrifugal Speed	220 rpm, 250 rpm, 300 rpm, 350 rpm
Electric Current	0.5A
Weight	1.2 kg
Voltage	220V
Power	50 Hz



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)