

T-MANIFOLD FREEZE DRYER LTFDR10-5



T-MANIFOLD FREEZE DRYER LTFDR10-5

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

LTFDR10-5 T-MANIFOLD FREEZE DRYER

7 Inch color touch control LCD display

Display run time & vacuum parameter

The cold trap with freezing function eliminates the low temperature refrigerator

Imported compressor, high capability of water-capture, low noise

Mixed refrigerant, green environmental protection

Built-in heat exchanger effectively reduces the temperature of cold trap

Anti-return oil vacuum pump to prevent the pollution of the material

Transparent bell-cap drying room is safe and visual

Optional distribution defrosting function

Optional nitrogen inflating valve



SPECIFICATIONS

Model	LTFDR10-5
Type	T-Manifold
Plate Load Capacity	1 L
Condenser Temperature	<50°C
Water Holding Capacity	4 Kg/24 h
Freeze Drying Surface Area	0.092 m ²
Vacuum Degree	<10 Pa
No of Sealed Valves	24 pieces for 24 ampoules
Overall Dimension (mm)	500x400x650 mm
Power	950 W
Power Supply	220 V 50 HZ /110 V 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com