

# TOP PRESS FREEZE DRYER LTFDR13-15



# TOP PRESS FREEZE DRYER LTFDR13-15

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## LTFDR13-15 TOP PRESS FREEZE DRYER



### SPECIFICATIONS

|                       |                                       |
|-----------------------|---------------------------------------|
| Model                 | LTFDR13-15                            |
| Type                  | Top Press type electric heating shelf |
| Condenser Temperature | -110°C                                |
| Freeze Drying area    | 0.08 m <sup>2</sup>                   |
| Bulk capacity         | 0.8 L                                 |
| Tray                  | φ180mm, 3pcs                          |



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)