

T-MANIFOLD FREEZE DRYER LTFDR13-26



T-MANIFOLD FREEZE DRYER LTFDR13-26

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

LTFDR13-26 T-MANIFOLD FREEZE DRYER



SPECIFICATIONS

Model	LTFDR13-26
Type	Manifold Type
Condenser Temperature	-110°C
Freeze Drying area	0.18 or 0.27 m ²
Bulk capacity	1.8 or 2.7 L
Tray	φ240mm, 4 or 6pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com