

# STANDARD FREEZE DRYER LTFDR13-31



# STANDARD FREEZE DRYER LTFDR13-31

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## LTFDR13-31 STANDARD FREEZE DRYER



### SPECIFICATIONS

Model	LTFDR13-31
Type	Standard type electric heating shelf
Condenser Temperature	-110°C
Freeze Drying area	0.18 or 0.27 m <sup>2</sup>
Bulk capacity	1.8 or 2.7 L
Tray	φ240mm, 4 or 6pcs



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)