

# TOP PRESS FREEZE DRYER LTFDR13-37



# TOP PRESS FREEZE DRYER LTFDR13-37

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## LTFDR13-37 TOP PRESS FREEZE DRYER



### SPECIFICATIONS

Model	LTFDR13-37
Type	Top Press with Manifold type electric heating shelf
Condenser Temperature	-110°C
Freeze Drying area	0.09 m <sup>2</sup>
Bulk capacity	0.9 L
Tray	φ200mm, 3pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)