

IN SITU FREEZE DRYER LTFDR14-5



IN SITU FREEZE DRYER LTFDR14-5

Also known as Silicone Oil Heating Lyophilizer.

LTFDR14-5 IN SITU FREEZE DRYER



SPECIFICATIONS

Model	LTFDR14-5
Type	Silicone Oil Heating In-situ pilot freeze Dryer
Condenser Temperature	-50°C
Drying area	0.2m ²
Tray temperature	-40°C to +50°C
Trays load capacity	2~3kg/24hr



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com