

# IN SITU FREEZE DRYER LTFDR14-8



# IN SITU FREEZE DRYER LTFDR14-8

Also known as Silicone Oil Heating Lyophilizer.

## LTFDR14-8 IN SITU FREEZE DRYER

One-piece design, compact equipment installation

The maintenance orientation is concentrated on one side, which is convenient for maintenance

A variety of security interlocking functions to fully protect the product

The equipment is conveniently connected to the clean room, and can be connected to the isolator

Realize dual control of clean room and maintenance room, but interlocking with each other

Simple operation, all run automatically



## SPECIFICATIONS

Model	LTFDR14-8
Type	Pilot In-situ freeze Dryer
Condenser Temperature	-75°C
Drying area	0.24m <sup>2</sup>
Trays load capacity	6.5kg/24hr
Tray temperature	-50°C to +60°C
Computer connection	RS232

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)