

TOP PRESS WITH MANIFOLD FREEZE DRYER LTFDR9-4



TOP PRESS WITH MANIFOLD FREEZE DRYER

LTFDR9-4

Also known as Benchtop Freeze Dryer, Tabletop Lyophilizer, Laboratory Benchtop Freeze Dryer, Laboratory Tabletop Lyophilizer, Benchtop Lyophilizer, Laboratory Benchtop Lyophilizer, Top Press with Manifold Freeze Dryer.

LTFDR9-4 TOP PRESS WITH MANIFOLD FREEZE DRYER

- 7 Inch color touch control LCD display
- Display run time & vacuum parameter
- The cold trap with freezing function eliminates the low temperature refrigerator
- Imported compressor, high capability of water-capture, low noise
- Mixed refrigerant, green environmental protection
- Built-in heat exchanger effectively reduces the temperature of cold trap
- Anti-return oil vacuum pump to prevent the pollution of the material
- Transparent bell-cap drying room is safe and visual
- Optional distribution defrosting function
- Optional nitrogen inflating valve



SPECIFICATIONS

Model	LTFDR9-4
Type	Top Press with Manifold
Drying area	0.092 m ²
Tray	200 x 3 mm
Tray capacity Condensor	1 L
No-load temperature	< -55°C (-80 optional)
Vacuum	< 10 Pa
Interlayer spacing	70 mm
Vials number for reference	165 (22) 285 (16) 560 (12)
Ability to capture Water	3 kg
Machine Size	500x400x370 mm
Power	850 W
Voltage	220V/50Hz or 110V/60HZ



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com