

T-MANIFOLD FREEZE DRYER LTFDR9-5



T-MANIFOLD FREEZE DRYER LTFDR9-5

Also known as Benchtop Freeze Dryer, Tabletop Lyophilizer, Laboratory Benchtop Freeze Dryer, Laboratory Tabletop Lyophilizer, Benchtop Lyophilizer, Laboratory Benchtop Lyophilizer, T-Manifold Freeze Dryer.

LTFDR9-5 T-MANIFOLD FREEZE DRYER

7 Inch color touch control LCD display

Display run time & vacuum parameter

The cold trap with freezing function eliminates the low temperature refrigerator

Imported compressor, high capability of water-capture, low noise

Mixed refrigerant, green environmental protection

Built-in heat exchanger effectively reduces the temperature of cold trap

Anti-return oil vacuum pump to prevent the pollution of the material

Transparent bell-cap drying room is safe and visual

Optional distribution defrosting function

Optional nitrogen inflating valve



SPECIFICATIONS

Model	LTFDR9-5
Type	T-Manifold
Plate Load Capacity	1 L
Condenser Temperature	<50°C
Water Holding Capacity	3 Kg/24 h
Freeze Drying Surface Area	0.092 m ²
Vacuum Degree	<10 Pa
No of Sealed Valves	24 pieces for 24 ampoules
Overall Dimension (mm)	500x400x370 mm
Power	850 W
Power Supply	220 V 50 HZ /110 V 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com