

-10°C TO -25°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE18-4



-10°C TO -25°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE18-4

LTFRE18-4 -10°C TO -25°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT

Double Door Drawer x 12 USB port



SPECIFICATIONS

Model	LTFRE18-4
Capacity	358 L
Temperature Range	-10°C to -25°C
Alarm	High&low temp, power failure, control error, door open, high ambient temperture
External Dimension	720×830×1870 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com