

-10°C TO -40°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE22-7



-10°C TO -40°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE22-7

LTFRE22-7 -10°C TO -40°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT

Double Door, Drawer x 12, USB port



SPECIFICATIONS

Model	LTFRE22-7
Capacity	358 L
Temperature Range	-20 to -40°C
Alarm	High&low temp, power failure, control error, door open, high ambient temperature
External Dimension	720×830×1870 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com