

-40°C TO -86°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE27-16



-40°C TO -86°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE27-16

LTFRE27-16 -40°C TO -86°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT

Self-cascade system with mixed-refrigerant



SPECIFICATIONS

| | |
|--------------------|--|
| Model | LTFRE27-16 |
| Capacity | 728 L |
| Temperature Range | -40°C to -86°C |
| 2ml Sample Loading | 50400pcs |
| Alarm | High&low temp, power failure, control error, door open, abnormal voltage, low battery, condensor clean, high ambient temperature |
| Alarm port | USB+Remote alarm port |
| External Dimension | 1229×1012×1994.5 mm |
| Optional Parts | 110V transformer; Datalogger; Cold storage rack, box; LN2/CO2 backup system |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com