

-40°C TO -86°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE27-17



-40°C TO -86°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT LTFRE27-17

LTFRE27-17 -40°C TO -86°C ULTRA LOW TEMPERATURE FREEZER UPRIGHT

Cascade system, HC green refrigerant R290 & R170, Energy saving), Touch Screen, VIP



SPECIFICATIONS

Model	LTFRE27-17
Capacity	728 L
Temperature Range	-40°C to -86°C
2ml Sample Loading	50400pcs
Alarm	High&low temp, power failure, control error, door open, abnormal voltage, low battery, condensor clean, high ambient temperature
Alarm port	USB+Remote alarm port
External Dimension	1229×1012×1994.5 mm
Optional Parts	110V transformer; Datalogger; Cold storage rack, box; LN2/CO2 backup system



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com