

PORTABLE FREEZER LTFRE34-2



PORTABLE FREEZER LTFRE34-2

LTFRE34-2 PORTABLE FREEZER

Real time temperature value display

Refrigeration and temperature control to ensure the consistency of samples

Easy to operate and could be used repeatedly

Meet the basic requirements of cold chain transportation and storage of vaccines



SPECIFICATIONS

| | |
|-----------------------|---|
| Model | LTFRE34-2 |
| Inner Capacity | 8 L |
| Outer size | 335x225x205 mm |
| Product Configuration | 1 box 1 digital thermometer 2 95KPa biosafety transport tanks 2 sponge reagent tubes 2 pieces of absorbent cotton 4 crystal ice boxes |
| Material | PP |
| Insulation Material | EPS |
| Gross weight | 2.5 kg |
| Packing size | 380x260x290 mm |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com